

# Bebe-Ball



- Bebe-Ball**
- Bulb suction adaptor for Bebeonkers 16-1 and 16-2
  - Helps prevent damage that occurs if one-hole device accidentally touches the sensitive mucosa of a newborn or preemie
  - Disposal after 12 to 24 hours or per hospital policy to prevent infection is recommended
  - Use on transports or when power is unavailable and STAT suctioning is required
  - Standard bulb syringe has one hole and is inappropriate because it can grab and tear tender tissue too


Bebe-Ball Manual Suction Device		
Order #	Product Description	Bx Qty
Bebe-Ball with Bebeonkers		
16-4	Bebe-Ball with Bebeonkers 16-2 (Std)	50
16-5	Bebe-Ball with UP Bebeonkers 16-1 (bulb)	50
Bebe-Ball with Bebeonkers and Keeper Cover		
16-4KK	Bebe-Ball with Bebeonkers 16-2 (Std) with Keeper Cover	50
16-5KK	Bebe-Ball with UP Bebeonkers 16-1 (bulb) with Keeper Cover	50



Red ball is our trial color, actual production will be blue translucent.



# Bebe Keeper Kover




Can be used up to 12 hours on the same patient, or per hospital policy.



Bebe Keeper Kover is an effective cover to keep the Bebeonkers clean.

Can be used on Bebeonkers I and Bebeonkers II with suction on or off, protecting the device while not in use.



Order #	Product Description	Bx Qty
<b>Bebe Keeper Kover Only</b>		
16-3	Bebe Keeper Kover for Bebeonkers 6 pcs/bag, 50 bgs/cs	300

Can also be used with the Bebe-Ball.



Note new blue translucent color