

Better for the Baby

Better for the Nurse

Better For The NICU/Hospital

Better X 3

CAT



















www.small-beginnings.com 1-800-676-0462 Fax orders to 1-760-948-1916



This page is dedicated to Nurses, who have supported us and helped us bring our products to the neonatal intensive care units around the world. Your desire to do your best for the infants in your care has been the reason Small Beginnings Inc. is considered the leader in producing great patented neonatal devices. When you called us with needs or problems about products that didn't work, we have rolled up our sleeves and created solutions.

Our Mission Statement is as true today as it was when we started our company.

We will make products that are better for the Baby! We will make products that are better for the Nurse! We will make products that are cost-effective and better for the hospitals! Better x 3.

It is this philosophy that has individual hospitals and entire hospital systems converting to our wonderful product line. We are proud to bring you the very best year after year and with your support we will continue to do so.

My Sincere Thanks

M. Sharon Rogone RNC-E Chief Executive Officer Small Beginnings Inc. A Certified Women



Table of Contents				
Category & Hem Pg # Category & Hem				
Fantastic Four		*continued	1000	
Bílí-Bonnet	4	Reflux Wedge	24	
Pacee Pacifier	30	Bed Buddy	26	
Cuddle Buns Diapers	42	Diana Bed Buddy	28	
Bebeonkers 0, 1, 11, 111	32	Other solutions		
Suctioning Solutions		Climate Cover	6	
Belleonkers 0, 1, 11, 111	32	Quillo-iso Cover	2	
Bebe-Ball	34	Pacee Pacifier	30	
Bebe-Keeper Kover	35	Parent Guides	36	
Positioning Solutions		Preemie Gowns	37	
Heavenly Foam Topper	10	Wing-slings	38	
Wonder Star	12	Privacy Poncho	39	
Cuddle-up	14	Dríp-drop Solutions		
Nests	16	Bebe-Chix	40	
Positioner Sampler	12	Birth N Bath Wraps/Bags	41	
Píllow Sampler	19	Cuddle Buns Diapers	42	
Posture Pillows lab, head wedgel	20	Beary Bili Buns Diapers	44	
Heck rolls	22	Notes	46	

Bili-Bonnet





Stretches open to allow Bili-light to penetrate the entire head.

3 Sizes Fit all head shapes.

- · No Confusion
- · Easy to apply
- · No straps to frustrate
- No adhesives to damage skin



- Compression Technology

 Makes ours the only mask with
 minimal light penetration
- Rests on ocular face bones and removes pressure from eye
- Mask is Reversible made of PE Foam & Brushed Nylon
- Conforms to baby's face with plenty of room for adjustment on each side

Order #	Product Description	Bx Qty
	Bili-Bonnet Phototherapy Masks	
01A	Preemie size for infants approx. 27-30 weeks	50
018	Newborn size for infants above 31 weeks approx.	50
01 <i>C</i>	Ultra preemie size for infants below 27 weeks approx.	50

Tabs are rounded & very low profile hook & loop that allow attachment at any angle to maximize fit and minimize slippage.

Bili-Bonnet

Latex Free-No adhesives

How to put the patented Bili-Bonnet phototherapy mask on correctly



Remove from package, then stretch gauze and pull over fist to create a bonnet.



Stretch bonnet over baby's head to go from the nape of the neck to the forehead.



Attach velcro tab to one side.

Note, the mask has a 1 inch
area for better adjustment &
velcro tabs are rounded.



Repeat with tab on other side. The mask should be snug & secure but not tight.



Ensure the netting is opened up all the way over the head to allow light to penetrate.



Now once again stretch netting down to the nape of the neck and up to the forehead.



Most Important Step
Once all is in place and
secure, unhook and reattach
the velcro tabs to pull
up/down/midline as head
shape dictates.



The end result is a mask that stays on better giving the nurse more time for other things.

Protects the eyes better & allows more light to the head.

Take time to save time!


Climate Cover III



- "The perfect micro-environment"
- Reduces insensible water loss and cold stress
- Reduces energy cost as it maintains body temperature
- Becomes an OXY-BODYHOOD when infant is off the ventilator
- **Warning must use with O2**

Order #	Product Description	8× Qty
95 AV	Climate Cover III	
02E	The imoregulation Oxy-Body Hood	6
	19×25×10	

- · Can stay in place for most care with easy access Arm Ports
- · Nurses can do good developmental care and positioning
- · Can be X-rayed through
- Maintains a neutral thermal environment even if/when heat source must be removed for a brief time





For removal in case of emergencies or other procedures all lines and tubing go under the bottom arches on all 4 sides

Climate Cover III





Can make every radiant warmer a

"Poor-man's Giraffe"



Iso-Sphere "Quillo"

Universal Isolette Cover

Infants in the NICU are bombarded with constant noise, harsh lights and painful procedures. These tiny babies learn quickly that when they are covered they do not get poked, prodded and stuck. Reducing stress promotes faster healing, ideal growth and proper development which in turn reduces hospital stays with better outcomes.

- 8 ounce guilting to reduce noise & light
- Machine washable and reusable
- Cotton/poly blend & flame resistant
- Gray non-stimulating inside reduces stress and promotes better rest/sleep These covers serve a great purpose!!





Order #	Product Description	B× Qty
	Iso-Sphere Cover	Î
	Universal Isolette Cover	

Fold into a pocket "Quillo"

These flame resistant covers not only fit over any isolette, bassinette or even cribs, they also easily cover the "Giraffe" as well. When not in use it folds into itself to make a pillow or "Quillo". It can also be labeled with your hospital and/or the baby's name and birth date for a perfect keepsake from the NICU.

Can also be carried in your hospital's gift shop.



Heavenly Foam





2 Straps to
Use as Needed

Recommend Use with Bed Buddy

Heavenly Foam
Included

Visco-Elastic, first developed by NASA

Latex free memory foam

Perfect topper for the neonatal bed

Reduces surface tension, pressure points & stress

In 2 sizes and is 5/8 of an inch thick



Heavenly Foam is included and used with Cuddle Up. Note use with other products as well as Bed Buddy to provide needed soft rebound and boundaries.




Vinyl Covered Heavenly Foam with Knit Cover

	Heavenly Foam Mattress	
Order #	Product Description	Cs Qty
	Heavenly Foam Mattress only	
04VE	Preemie/Newborn Heavenly Foom Mattress - 20"x13.5"x5/8"	24
04CVE	Ultra Preemie Heavenly Foam Mattress 10"×13.5"×5/8"	10
	Leatherette Vinyl Covered Heavenly Foam	
04VE-C	Preemie/Newborn Vinyl Covered mattress	В
04CVE-C	Ultra Preemie Vinyl Covered mattress	6
04A VE-C	Vinyl Covered mattress for Preemie/Newborn with Fleece Covers	8
04AKVE-C	Vinyl Covered mattress for Preemie/Newborn with Knit Covers	В
04CCVE-C	Vinyl Covered mattress for Ultra Preemie with Fleece Covers	6
04CCKVE-C	Vinyl Covered mattress for Ultra Preemie with Knit Covers	6

"We only make products that are better for the baby, better for the nurse and better for the hospital." – M. Sharon Rogone RNC-E CEO

Heavenly Foam Covers



- · Reece, knit & disposable
- Waterproof disposable
- Vinyl covered Heavenly Foam
 also available





3	Waterproof	Disposab	le Cover
	200		1
		1	1
		14.236	1
	The same of		
160			11
	and the same of th		

Order #	Product Description	Cs Qty
	Covers for Vinyl Covered Heavenly Foam Mattress	
04A	Preemie/Newborn Fleece Covers	10
04AK	Preemie/Newborn Knit Covers	10
04C	Ultra Preemie Pleece Covers	10
0.4CK	Ultra Preemie Knit Covers	10
04VED-C	Preemie/Newborn Disposable Covers	20/Pkg
04CVED	Ultra Preemie Disposable Covers	20/Pkg
	Disposable Covers for Heavenly Foam Mattress	
04VED	Disposoble covers for Preemie/Newborn Mattress	20/Pkg
04CVED	Disposable covers for Ultra Preemie Mattress	20/Pkg
	Baby Blankee	
04b	Boby Blonket 25"x31"	12

Wonder Star

Wonder Star

- Multi-purpose positioning device
- Soft gel rebound beads
- Washable
- Move & manipulate in many ways
- Keep baby comfortable.
- 2 sizes
- Fleece or Knit
- Great addition to any NICU.



Order #	Product Description	Cs Qty
Wonde	er Star Multi-Purpose Positioning Pillo	ow/
05A	Wonder Star in Fleece 15"arm span and 8" wide body	10
05AK	Wonder Star in Knit 15"arm span and 6" wide body	10
050	Lil' Wonder Star in Fleece 13"arm span and 6.5" wide body	10
05 <i>c</i> k	Lil' Wonder Star in Knit 13"arm span and 6.5" wide body	10



Wonder Star



Very flexible positioning

Available in Fleece or Knit this multi-functional positioning aid can meet many needs. Providing the best positioning for your baby, whether maintaining head or body position or nestling baby in a womb-like position.

The Wonder Star can do it all.





Cuddle Up

The Wonderful Cuddle Up

Comes with a Heavenly Foam

mattress to put in its

pillow case sleeve.

Recommend use with the

Bed Buddy to provide soft,

push back boundaries

**Ultra Preemie ames with Bed Buddy

Nurses requested a better swaddling device, so Sharon Rogone created the Cuddle Up Swaddler.

		THE PERSON NAMED IN
Order #	Product Description	Cs Qty
	Cuddle Up Sets	
D6A	Preemie/Newborn Fleece Cuddle Up with Foam Mattress	6
D6AK	Preemie/Newborn Knit Guddle Up with Foam Mattress	б
060	Ut ra Preemie Cudd te Up with Faam Mattress and 30°Bed Buddy	LO
D6CK	Uttra Preemie Cudd le Up writh Faam Mattress and 30° Bed Buddy	LO
	Cuddle Up Snugglers	
065 A	Preeme Cuddle Up Snuggler	6
065 B	Newbarn Cuddle Up Snuggler	6
065 C	Ultra Preemie Cuddle Up Snuggler	6
065 CW	Ultra Preemie Cudd le Up Snuggle r Warming Winap	6

Introducing the new Snuggler wrap
designed to maintain infants
temperature and secure tubing and
lines during any procedure; such as,
transports, weighing, exams, and
kangaroo care, etc. Is also available
in ultra preemie size only, with
2 pockets on the back side for
placement of heel warmers.





Cuddle Up Covers

Fleece Cuddle Up Cover



Note 2 straps

With inner sleeve for the Heavenly Foam Mattress.



Disposable Cuddle Up Cover


Order #	Product Description	Ćs Qty
	Cuddle Up Covers only	
06AA	Preemie/Newborn Fleece	10
06AAK	Preemie/Newborn Knit	10
06CC	Ultra Preemie Fleece	10
06CCK	Ultra Preemie Knit	10
06AAD	Preemie/Newborn Disposable Cover	20/Pkg
0.6CCD	Ultra Preemie Disposable Cover	20/Pkg

Cuddle Up Covers

- Fleece, Knit or Disposable Covers
- 2 straps to Cuddle, snuggle and secure
- Sleeve for Heavenly Foam Mattress

How do you keep the Heavenly Foam mattress pad from being soiled or wet?
The answer is our waterproof disposable covers lined with plastic to fix the problem.
Now your baby can lay on NASA's best for comfort and circulation
and know it is clean and dry.

Waterproof Disposable Cover





Nest Sets

Created for the sensitive, hard to manage infant

- Decrease decels and desats
- Reduce stress and negative reactions & set backs
- Gel beads push back for spatial & brain development
- Bags support baby's back, front or side
- · Provides for proper positioning
- 3 sizes: Preemie, Newborn and Ultra Preemie to meet the needs of all age & size babies



Super soft gel beads

Patented by our founder, the Preemic Nest replaces rolled blankets & time lost by the nurse



The ultimate swaddler

We added a 4th sleeve so you can rearrange the bags & only turn baby

Now with Velcro tabs to keep bags from sliding out



Order #	Product Description	Cs Qty
	Nest Sets Infant Positioning Device	
09C	Preemie Nest Set with Fleece Pad	10
09CK	Preemie Nest Set with Knit Pad	10
09CC	Newborn Nest Set with Fleece Pad	6
09CCK	Newborn Nest Set with Knit Pad	6
09CCC	Ultra Preemie Nest Set with Fleece Pad	10
09CCCK	Ultra Preemie Nest Set with Knit Pad	10

How to prepare the Nest

- Center baby in supine position
- Next grab bead bag close to you
- Roll the bag toward yourself to turn baby from supine to side position
- · Push other bead bags in against baby
- Surrounds baby with womb-like feeling, soft, cushy, rebound

Turning & Positioning in seconds!!










The Nest Pads are designed with ribbed knit sleeves sewn into the back of a super soft fleece or knit blanket. Each sleeve then has Velcro at each end to keep the bags from falling out. Pads come with two long and two short sleeves for ease in turning and positioning. Disposable pads are also available



Velcro fastening to hold bead bags in place

Order #	Product Description	Cs Qty
	Nest Pads only	5
09A	Preemie Pleece Pads	12
09AK	Preemie Knit Pads	12
09AA	Newborn Fleece Pods	12
09AAK	Newborn Knit Pads	12
09AAA	Ultra Preemie Fleece Pads	12
09AAAK	Ultra Preemie Knit Pads	12
	Disposable Pads Only	Pkg Qty
09AAD	Preemie/Newborn Disposable Pads	20
09AAAD	Ult ro. Preemie Disposable Pads	20





All the positioning aids
any unit could ever want or need.
Every situation met and enough
for more than one baby.
All bags wipe clean, are reusable
and come with washable or
disposable covers.
Guaranteed for 2 years.
Covers come in knit, fleece or
disposable material.

All the amazing benefits of our 7 positioning products.
Assists PT, OT & all, with positioning, alignment, comforting and swaddling.
Remember when we position, comfort and reduce stress we are enabling normal growth and development.
When we mimic the womb emironment we help the cognitive, spatial and physical growth of our tiny charges.
Never forget our oath to "Do No Harm."

Order #	Product Description	Cs Qty
	Positioning Sampler Kits	
07A	Preemie Positioning Sampler Kit includes 1 of each device with fleece covers: 10A Abdominal Pillow, 10B Head Pillow, 10C Wedge Pillow, 09C Nest Set,06A Cuddle-Up, 05A Wonder Star, 14A Bed Buddy.	1 each/ 1 set
07AK	Preemie Positioning Sampler Kit includes 1 of each device with knit covers: 10A Abdominal Pillow, 10B Head Pillow, 10C Wedge Pillow, 09C Nest Set, 06A Cuddle-Up, 05A Wonder Star, 14A Bed Buddy.	1 each/ 1 set
07C	Ultra Preemie Positioning Sampler Kit with 1 of each device with fleece covers: 10AA Abdominal Pillow, 10BB Head Pillow, 10CC Wedge Pillow, 09CCC Nest Set, 06C Cuddle-Up, 05C Wonder Star, 14E Bed Buddy.	1 each/ 1 set
07CK	Ultra Preemie Positioning Sampler Kit with 1 of each device with knit covers: 10AA Abdominal Pillow, 10BB Head Pillow, 10CC Wedge Pillow, 09CCC Nest Set, 06C Cuddle-Up, 05C Wonder Star, 14E Bed Buddy.	1 each/ 1 set

Posture Pillow Sampler Kits

***Reference Article
Premature Infants: A Later Look
NYU Child Study Center


"Positioning the preemie can also influence future development. For example, preemies may be routinely left to lie on their backs, a position that differs from the natural curved fetal position they would maintain in utero. This can result in a shortening of the scapula muscles between the shoulder blades, causing a restricted range of motion and later motor problems. Thus, the older infant may be unable to appropriately reach for objects in front of him..."



Posture Pillow Sets Includes 1 of each pillow:

- Abdominal/prone positioning
- Head/neck or head
- Wedge/all around
- With Knit or Pleece covers









Disposable





2	C .	4	1 1 1
-	Sizes	ANTI	മവഭ

Order #	Product Description	Cs Qty
	Posture Pillow Sets with Gel Beads	
10b	Set: 1 Preemie Ab, Head, & Wedge Pillow & Pleece Covers	10 Sets
10DK	Set: 1 Preemie Ab, Head, & Wedge Pillow & Knit Covers	10 Sets
1000	Set: 1 Ultra Preemie Ab Head & Wedge Pillow & Fleece Covers	10 Sets
10bbk	Set: 1 Ultra Preemie Ab, Head & Wedge Pillow & Knit Covers	10 Sets

Posture Pillows

Gal Baseds or Visco Bastic (VE)

Order #	Product Description	Cs Qty
t t	Foam Posture Pillows with Covers	
10A-VE	Preemie Abdominal Pillow with Fitted Fleece Cover	10
10AK-VE	Preemie Abdominal Pillow with Fitted Knit Cover	10
10AA-VE	Ultra Preemie Abdominal Pillow with Fitted Fleece Cover	10
10AAK-VE	Ultra Preemie Abdominal Pillow with Fitted Knit Cover	10
108-VE	Preemie Head Pillow with Fleece Cover	10
108K-VE	Preemie Head Pillowwith Knit Cover	10
1088-VE	Ultra Preemie Head Pillow with Fleece Cover	10
1088K-VE	Ultra Preemie Head Pillowwith Knit Cover	10
108NB-VE	Newborn Head Pillow with Fleece Cover	10
108NBK-VE	Newborn Head Pillow with Knit Cover	10



X-Ray thru

Gel Bead Pillows



Head Pillow-Gel Beads or VE Foam

- Allow better shaping of head
- Beads can simply be moved to one side to make a neck roll
- Thick VE memory foam you can x-ray through



Abdominal Pillow-Gel Beads or VE Foam

- Softer TPE beads that feel like gel but never evaporate
- · High profile VE foam for greater height and can be x-rayed thru
- · Both developmentally designed with fitted covers
- Better prone positioning, better lung expansion & oxygenation

	Boad Bag Pasturo Pillows with Covers	
LDA	Preemie Abdamma Pillaw with Fitted Fleece Caver	LD
LOAK	Preemie Abdammal Pillaw with Frited Knit Caver	LD
LDAA	UhraPreemie AbdomnalPillow with Fitted Fleece Cover	10
LDAAK .	Ultra Preeme Abdomnal Pillow with Frited Knit Cover	LO
108	Preemie Head Pillow with Fleece Cover	LD
LOOK	Preemie Head Pillow with Knit Caver	LD
1088	Utira Preemie Head Pillow with Fleece Cover	LD
LOBOK	Ultra Preemie Head Pillow with Knit Cover	10
100	Preemie Wedge Pillow with Fleece Cover	10
LOCK	Preemie Wedge Pillow with Knit Cover	LO
LDCC	Ultra Preemie Wedge Pillow with Fleece Caver	LD
LDCCK	Uhra Preemie Wedge Pillaw with Knit Caver	LD



Wedge Pillows

- Multi-Functional
- Holds vent tubing & IV lines
- Quickly changes baby's position
- Wedge Pillow is a work horse
 You're sure to find your own uses!

Posture Pillow Covers

Fleece **

- Very soft
- Light-weight
- 100% polyester
- Mach ine washable
- Cheery colorful
- Multi-animal print

Knit **

- Very soft
- 50/50 cotton/poly
- Machine washable
- Cheery colorful
- Multi-animal print
- White background

Disposable **

- Soft
- Sturdy
- Non-woven
- Water resistant
- · Cheery colorful
- Multi-an imal print

Order #	Product Description	Pkg Qty
	Posture Pillow Covers Only	
10E	Preemie Head & Wedge Pillow Pleece Covers	10
10EK	Preemie Head & Wedge Pillow Knit Covers	10
10EA	Preemie Abdominal Pillow Pitted Pleece Covers	10
10EAK	Preemie Abdominal Pillow Pitted Knit Covers	10
10EE	Ultra Preemie Head & Wedge Pillow Fleece Covers	10
10EEK	Ultra Preemie Head & Wedge Pillow Knit Covers	10
10EEA	Ultra Preemie Abdominal Pillow Fitted Fleece Covers	10
10EEAK	Ultra Preemie Abdominal Pillow Fitted Knit Covers	10



- Ab cover is fitted to new updated shape
- Wedge coverdesigned to fit wedge pillow
- Head cover designed to fit head pillow



Order #	Product Description	Pkg Qty			
	Disposable Covers				
10ED	Preemie Head & Wedge Pillow Disposable Covers	50			
10EAD	Preemie Abdominal Pillow Disposable Covers	50			
10EED	Ultra Preemie Head & Wedge Pillow Disposable Cove	20			
10EEAD	Ultra Preemie Abdominal Pillow Disposable Covers	50			



Packaged Posture Pillow Disposable Covers

Neck Rolls

Undershirts no more!!!

Our soft gel-like beads

are flexible & conform

to fit comfortably & allow
neck to remain in a slightly
extended position for

easier ventilation.









Order #	Product Description	Cs Qty
	Bead Neck Rolls Infant Positioning Device	
10N	Presmis Neck Roll with Fleece Cover	10
10NK	Preemie Neck Roll with Knit Cover	10
10NLT	Late Term Preemie Neck Roll with Fleece Cover	10
10NLTK	Late Term Preemie Neck Roll with Knit Cover	10
10NN	Ultra Preemie Neck Roll with Fleece Cover	10
10NNK	Ultra Preemie Neck Roll with Knit Cover	10
10NB	Newborn Neck Roll with Fleece Cover	10
10NBK	Newborn Neck Roll with Knit Cover	10
	Foam Neck Rolls Infant Positioning Device	
10N-VE	Presmis Neck Roll with Fleece Cover	10
10NK-VE	Preemie Neck Roll with Knit Cover	10
10NLT-VE	Late Term Preemie Neck Roll with Fleece Cover	10
10NLTK-VE	Late Term Preemie Neck Roll with Knit Cover	10
10NN-VE	Ultra Preemie Neck Roll with Fleece Cover	10
10NNK-VE	Ultra Preemie Neck Roll with Knit Cover	10
10NB-VE	Newborn Neck Roll with Fleece Cover	10
10NBK-VE	Newborn Neck Roll with Knit Cover	10

Neck Roll Covers





Neck Rolls with vinyl wipe clean surface & very soft gel-like beads take the place of those rolled-up undershirts, etc. Available in 4 sizes to meet the needs of large to micro-size babies. A fleece cover is included, knit & disposable covers are also available. The knit & fleece covers are machine washable. Place Neck Roll into a cover and place behind infants neck in supine position, be careful not to hyper extend or impair circulation in the neck. These can be used as hand rolls, as well.



Disposabla covar

Fleece Cover



Knit Cover

Order #	Product Description	Cs Qty
	Neck Rolls Covers Only	
10NE	Preemie Neck Roll Fleece Cover	10
10NEK	Preemie Neck Rolls Knit Covers	10
10NLTE	Late Term Preemie Neck Rolls Fleece Covers	10
10NLTEK	Late Term Preemie Neck Rolls Knit Covers	10
10NNE	Ultra Preemie Neck Rolls Fleece Covers	10
10NNEK	Ultra Preemie Neck Rolls Knit Covers	10
10NBE	Newborn Neck Rolls Fleece Covers	10
10NBEK	Newborn Neck Rolls Knit Covers	10
	Disposable Covers Only	Pkg Qty
10NED	Preemie Neck Rolls Disposoble Covers	20
10NLTED	Late Term Preemie Neck Rolls Disposable Covers	20
10NNED	Ultra Preemie Neck Rolls Disposable Covers	20
10NBED	Newborn Neck Rolls Disposoble Covers	20

Prints may vary.

Reflux Wedge

The Reflux Wedge is made from our durable, wipe clean vinyl with comfortable VE foam.
One side is for ultra preemies
And the other for preemies.
Fits in most isolettes or bassinets.

With Preemie size infant and Disposable Cover





	Bed Buddy Infant Positioning Devices	
Order #	Product Description	Cs Qty
	Bed Buddy with Covers	
14A	Preemie 44" Bag with Fleece cover	6
14AK	Presmis 44" Bog with Knit cover	6
148	Newborn 60" Bag with Fleece cover	6
148K	Newborn 60" Bag with Knit cover	6
14C	Ult ra Preemie 36" Bag with Fleece cover	10
14CK	UltraPreemie 36"Bagwith Knit cover	10
14E	XS Ultra Preemie 30" Bag with Fleece cover	10
14EK	X5 Ultra Preemie 30" Bag with Knit cover	10
	Bed Buddy Bag only	
14A-B	Preemis 44"Bag	6
148-8	Newborn 60"Bag	6
14C-B	Ultra Preemie 36" Bag	10
14E-B	XS Ultra Preemie 30" Bag	10

Reflux wedge Covers

The Diana Bed Buddy covers come in 3 very soft fabrics, fleece, knit and disposable. The fleece and knit can be machine washed and dried.







	Bed Buddy Infant Positioning Devices	
Order #	Product Description	Cs Qh
	Bed Buddy with Covers	
14A	Presmis 44" Bag with Fleece cover	6
14AK	Preemie 44" Bag with Knit cover	6
148	Newborn 60" Bag with Fleece cover	6
148K	Newborn 60" Bag with Knit cover	6
14C	Ultra Preemie 36" Bag with Fleece cover	10
14CK	Ultra Preemie 36" Bag with Knit cover	10
14E	X5 Ultra Preemie 30"Bag with Pleece cover	10
14EK	X5 Ultra Preemie 30" Bag with Knit cover	10
- "	Bed Buddy Bag only	
14A-8	Preemie 44"Bag	6
148-8	Newborn 60"Bag	6
14C-B	Ultra.Preemie 36"Bag	10
14E-8	XS Ultra Preemie 30" Bag	10

Bed Buddy

The Bed Buddy is the easiest positioning pillow to use. It quickly surrounds and interfaces with infant, providing a soft rebounding surface. It has zero replacement cost. Just wipe vinyl clean with any disinfectant, wash and replace cover and it's ready to go.






Wipe clean vinyl bead bag

	Order #	
1 S 1 S 2		
	144	
	14AK	
	148	
I. Marin Mar	148K	
	14C	
	14CK	
7	14E	X
	14EK	3
50		
	144-8	
	148-8	
	14C-B	
	14E-8	

	Bed Buddy Infant Positioning Devices	
Order #	Product Description	Cs Qty
	Bed Buddy with Covers	
144	Preemie 44" Bag with Fleece cover	6
14AK	Preemis 44" Bog with Knit cover	6
148	Newborn 60" Bag with Fleece cover	6
148K	Newborn 60" Bag with Knit cover	6
140	Ultra Preemie 36" Bag with Fleece cover	10
14CK	Ultra Preemie 36" Bag with Knit cover	10
14E	XS Ultra Preemie 30"Bag with Fleece cover	10
14EK	XS Ultra Preemie 30" Bag with Knit cover	10
	Bed Buddy Bag only	
14A-B	Preemie 44"Bag	6
148-8	Newborn 60"Bag	6
14C-B	Ultra Preemie 36" Bag	10
14E-8	XS Ultra Preemie 30" Bag	10

Bed Buddy Covers





Fleece Cover

Order #	Product Description	Cs Qty
	Bed Buddy covers only	
1444	Presmis Flasce Covers	10
14AAK	Preemie Knit Covers	10
1488	Newborn Fleece Covers	10
1488K	Newborn Knit Covers	10
1400	Ultra Preemie Fleece Covers	10
14CCK	Ultra Preemie Knit Covers	10
14EE	XS Ultra Preemie Fleece Covers	10
14EEK	XS Ultra Preemie Knit Covers	10
	Disposable Covers only	Pkg Qty
144 AD	Preemie Disposable Covers	20
1488b	Newborn Disposable Covers	20
14CCD	Ultra Preemie Disposable Covers	50
14EED	XS Ultra Preemie Disposable Covers	20




The Bed Buddy can be covered with our washable Fleece or Knit and is also available in our print disposable.

> All come with convenient pocket closure

Disposable Cover

Diana Bed Buddy

The Bed Buddy is the easiest positioning pillow to use.
It quickly surrounds and interfaces with infant,
providing a soft rebounding surface. It lasts and
has a 2 year guarantee. Just wipe vinyl clean
with any disinfectant, wash and replace cover
and it's ready to go.



	Bed Buddy Infant Positioning Devices	
inder #	Product Description	Cs Qt
	Bed Buddy with Covers	
14A	Preemie 44" Bag with Pleece cover	6
14AK	Preemie 44"Bag with Knit cover	6
148	Newborn 60" Bag with Fleece cover	6
148K	Newborn 60" Bag with Knit cover	6
14C	Ultra Preemie 36" Bag with Pleece cover	10
14CK	UltraPreemie 36"Bag with Knit cover	10
14E	XS Ultra Preemie 30"Bag with Fleece cover	10
14EK	X5 Ultra Preemie 30"Bag with Knit cover	10
	Bed Buddy Bag only	
14A-B	Preemie 44"Bag	6
14B-B	Newborn 60"Bag	6
14C-B	Ultra.Preemie 36"Bag	10
14E-8	XS Ultra Preemie 30"Bag	10

The Diana



Diana Bed Buddy Covers

The Diana with fleece cover.

The Diana Bed Buddy covers come in 3 very soft fabrics, fleece, knit and disposable. The fleece and knit can be machine washed and dried.







Note Cuffed End



	Bed Buddy Infant Positioning Devices			
Order #	Order # Product Description			
	Bed Buddy with Covers			
144	Preemie 44" Bag with Pleece cover	6		
14AK	Preemie 44"Bag with Knit cover	6		
148	Newborn 60" Bag with Fleece cover	6		
148K	Newborn 60" Bag with Knit cover	6		
14C	Ultra Preemie 36" Bag with Pleece cover	10		
14CK	UltraPreemie 36"Bag with Knit cover	10		
14E	XS Ultra Preemie 30"Bag with Fleece cover	10		
14EK	X5 Ultra Preemie 30"Bag with Knit cover	10		
	Bed Buddy Bag only			
14A-8	Preemie 44"Bag	6		
148-8	Newborn 60"Bag	6		
14C-B	Ultra.Preemie 36"Bag	10		
14E-8	XS Ultra Preemie 30" Bag	10		

Pacee Pacifier



- 100% Silicone all one piece as recommended by The Journal of Pediatrics.
- 6 Large ventilation holes makes it the safest pacifier at discharge and for home.
- 2 sizes in easy to see colors
- Teal for newborns and purple for preemies.






Order #	Product Description	B× Qty
	Pacee Pacifiers	
15A	Lil' Paces presmis pacifier 100% silicons	100
15B	Nu'Pacee newborn pacifier 100% silicone	100

Pacee Pacifier







Use shallow dip for CPAP

- Never needs cutting
- · Use deep V for ET Tube
- Use shallow dip for Nasal CPAP
- Collapses easily to the base
- · Mimics nutritive suck



- Softer nipple for premature infants
- · Restful non-nutritive suck
- · Won't easily fatigue infant
- Recommended by Lactation Consultants



Bebeonkers

Bebeonkers II (16-2) Our original patented design is an Oral/Nasal suction device that does not grab tissue and effectively removes the most tenacious secretions, quickly and without compromising the infant.

Now available: 16-0 Webeonker's for the ultra preemie infant, just like the 16-2 but smaller.

Packaged six per bag.





Flexible & Soft

Bebeonkers I (16-1) was designed for the Ultra Preemie infant to facilitate quick removal of secretions both nasally and orally but still retains a three-hole configuration for safe oral care without hurting baby.

Remember "Do No Harm"

Order #	Product Description	8x Qty
	Bebeenkers Onal/Nasal Suction Device Only	
16 1	Bebeankers [Uttra.Preeme Bubed Tp 6 pcs/z p bag,50 bgs/bx	300
16 2	Bebeankers [[Standard Tip 6 pcs/zip bag, 50 bgs/b×	300
16-0	Webcankers Utira Preemie 57D 1 p 6 pcs/xp bag. 50 bags/b×	300
16-1KK	16-1KK Bebeankers C with kover 3 Bebeankers + 3 cover/bag, 50 bgs/b×	
16-2KK	Bebeankers [[with kover 3 Bebeankers + 3 cove Wbag, 50 bgs/bx	300
16-0KK	Webeankers with kover 3 Bebeankers + 3 cove Wbag, 50 bgs/b×	300



Bebeonkers





Devices altered to facilitate the removal of oral secretions or designed ineffectively for tenacious secretions.

There decrees now shared the legisla aread to exactly as Name sale whatever they as to rector now effective.

Major liability

A ill above are single hole device will gob fissue when successed and that's ball.

* remember a single hole device will gob fissue when successed and that's ball.

Patented
2 or more holes
Bebeankers has 3
Does not grab tissue

No trauma No harm Because there has never been an appropriate device for suctioning intubated or newborn infants, doctors, nurses and therapists have been using and creating inferior one-hole devices. They endanger the infants and make the hospital and medical professional liable if an infant is injured.

The patented
"BEBEONKERS" is that right device!!
It is the right tool for the right job.
Two or more holes and a tape red
inner wall does the job quickly and
without tearing tissue.

"Sometimes you just can't stand by when things are being done that can hurt a baby. We saw a problem that was international and made a product to solve it!"

—Philip N. Rogone PA. RCP.

Bebe-Ball



Bebe-Ball Manual Suction Device				
Order #	Order # Product Description			
	Bebe-Ball with Bebeonkers			
16-4	Bebe-Ball with Bebeonkers 16-2 (Std)	50		
16-5	Bebe-Ball with UP Bebeonkers 16-1 (bub)	50		
16-6	Bebe-Ball with Webeonkers 16-0 (UP)	50		
B	ebe-Ball with Bebeonkers and Keeper Kove	r		
16-4KK	Bebe-Ball with Bebeonkers 16-2 (Std) with Keeper Kover	50		
16-5KK	Bebe-Ball with UP Bebeonkers 16-1 (bub) with Keeper Kover	50		
16- 6KK	Bebe-Ballwith Webeonkers 16-0 (UP) with Keeper Kover	50		



Bebe-Ball


- Bulb suction adaptor for Bebeonkers 16-0, 16-1 and 16-2
- Helps prevent damage that occurs if one-hole device accidently touches the sensitive mucosa of a newborn or preemie
- Disposal after 12 to 24 hours or per hospital policy to prevent infection is recommended
- Use on transports or when power is unavailable and STAT suctioning is required
- Standard bulb syringe has one hole and is inappropriate because it can grab and tear tender tissue too



Red ball is our trial color, actual production will be blue translucent.

Bebe Keeper Kover


Can be used up to 12 hours on the same patient, or per hospital policy.





Bebe Keeper Kover is an effective cover to keep the Bebeonkers clean.

Can be used on Bebeonkers I and Bebeonkers II with suction on or off, protecting the device while not in use.



Can also be used with the Bebe-Ball.

Note new blue translucent color

Parent Guides

- Written for parents in simple and easy to understand language
- NIDCAP inspired book developed by nurses
- Knowledge of baby's reactions, positive & negative to the NICU
- Empowers parents to be more involved
- Assists with discharge teaching for JCAHO requirements
- Parents are signed off as they learn how to care for their baby



Order #	Product Description	Bx Qty
	Parent Guide	
17 A	English Guide for Parents with infants in the NICU	50
17AW	English Guide for Parents with infants in the NICU with Hooks	50
178	Spanish Guide for Parents with infants in the NICU	50
178W	Spanish Guide for Parents with infants in the NICU with Hooks	50
170	French Guide for Parents with infants in the NICU	50
17CW	French Guide for Parents with infants in the NICU with Hooks	5 0



French Available Upon Request



All the best, SECU Staff in health freguessign

Now it's time to . . .

Have supplies ready bilegers, shather, buttles, etc.) Get a tay test and adjust it and know how to use it Pulcous a chance for my well haloy care and he that MD know her the staff income if you waste your bully loop commenced.

		_	
		-	
-			

Precious Preemie/Newborn Gowns

- The Precious Preemie/Newborn Gown is Perfect
- One piece gown with 3 Velcro openings to minimize handling
- Opens in the front or back
- Allow easy access to IV lines & electrode wires
- Machine washable cotton/poly blend, flame resistant



"As clinicians on a mission, seeing nurses cutting products to make them work is all the motivation we need to create something better."

-Kønneth Croteau RRT COO

Wing Sling



Remember those big babies born with a broken clavicle? Keep them safe and comfortable with this easy alternative to under shirts and safety pins.

- Soft padded brushed nylon & foam for the back of baby's neck
- Flexible stretch netting to create a sling to protect and cradle the affected arm
- Allows for natural movement while holding in a normal secure position







Order #	Product Description	Pkg Qty
	Wing 5 ling	
11A	Arm Immobilizer for broken clavicles	6

Privacy Poncho



Created for busy units

- Eliminates screens and their problems
- Provides privacy for Skin2Skin/Kangaroo care
- · Soft fleece material keeps baby warm
- 100% polyester, machine washable, flame resistant
- Large neck opening so nurse can evaluate infant positioning, lines and mom's breast feeding
- Strong Velcro attachments on both sides to hold IV lines and ventilator tubing
- Designed to fit any Mom or Dad

Order #	Product Description	Cs Qty
	Privacy Poncho	
12A	Privacy cover for breast feeding in the NICU	10
12AO	Privacy cover for breast feeding at home (available online)	10



BeBe Chix



Under Pads for the Wee Ones



These under pads are designed for the little messes from feedings, IV, blood draw drips or just easy clean up. Can be used to protect Heavenly Foam, Cuddle Up, Nest Sets and more. Save on laundering. Inexpensive and very handy.

Onder #	Product Description	Pkg Q ty		
	BeBe Chix Under Pads			
19A	Infant UnderPads 12"×10"	20		



Birth n' Bath Wraps

Hooded and/or Bubbled

Disposable insulated wrap for wet infants

- Use at birth and/or at bath time
- Warm & Insulate against insensible water lass and drafts
- To insulate, protect & swaddle infant.
- Safe, minimizing sticking to wet skin with bubble or soft non woven inside.
- Either Quilted, soft, fluffy, absorbent & insulated inner liner or inner side bubbled & insulated.
- Plastic waterproof outer layer.
- Disposable and economical.



Use

- Open wrap on warmer.
- Place infant in center.
- If hooded adjust to fit head.
- Pull net bandage over head to close and secure.
- Open just at umbilicus for line placement.

Order #	Product Description	Pkg Qty
	Birth n' Bath Wraps	
0.68nB	Preemie/Newborn Absorbent & Insulating Wraps at birth n bath	10
06BnB Hood	Preemie/Newborn Hooded Absorbent & Insulat- ing Wraps for use at birth n bath	10
06BnB UP	Ultra Preemie Hooded Absorbent & Insulating Wraps for use at birth n bath	10
0.68nB CLW	Ultra Preemie Clear Insulating Wraps at birth n bath	5

Cuddle Buns

What makes our diapers so special?



The answer

The narrow crotch.

Our pads are figure eight in shape.

This helps keep baby positioned to

prevent abduction.

**Ref: article: NYU Child Study Center Premature Infants: A Later Look

"C orrectly sized diapers is one such simple invention. For example, it was found that preemies who are put in diapers that are too large must keep their legs spread apart in order to accommodate the extra fabric. This resulted in incorrect turn out of the hips and thus later motor problems."





Order #	Product Description	Bx Qty
	Cuddle Buns Developmental Diapers	
21A	Preemie Size fit infants 2-4 lbs or 1300-2000 gms	360
21B	Large Preemie fit infants 4-6 lbs or 1850-2750 gms	360
21C	Ultra preemie fit infants under 2 lbs or 1300 gms	360

Cuddle Buns



- A crotch small enough to prevent abduction
- Double elastic that really works to prevent leaks
- A one way liner to keep baby dry





Note: Diaper is too big.

Legs are spread wide.

- No gel to swell
- You can easily get a urine specimen
- A cut out to keep diaper away from umbilical lines



"The making of a preemie diaper was very difficult but the dangers to the baby with legs abducted and plastic tearing their sensitive skin kept me on course to find a solution. It was a 5 year struggle to produce a better diaper."


—M. Sharon Rogone RN

Beary Small

Bili-Buns Phototherapy Diapers


The Beary Small Economical Diapers 22A-22B-22C are designed with a narrow crotch and made with SAP for those situations where extra absorbency is required. Made with hook and loop attachments and newly improved all white opaque outer layer and no tape to stick to sensitive skin.

And now our Beary Small Diapers are permeable by photo therapy lights.





Order #	Product Description	8x Qty
	Beary Small Economical Photo Therapy Diapers	
22A	Preemie fit infants 1.5-2.5 bs or 800-1300 gms	360
228	Lg. Preemie fit infants 2.5-4 lbs or 1300-1600 gms	360
220	Ultra preemie fit infants under 1.5 lbs or 800 gms	360







Beary Small

Bili-Buns Phototherapy Diapers

Beary Small Bili-Buns

- Still the same but now better
- Bili-Light penetrates diaper
- Now white opaque outer layer
- Hook & loop closure & no tape
- Economical
- Very absorbent with gel
- 3 preemie sizes
- Narrow crotch

Test showing light penetrating the Bili-Buns Diaper through the fluff.







Notes



Knowledge is Power

Contains a wealth of knowledge concerning developmental care and outcomes, also safe practices and standards of care for nurses; as well as the cost savings for the NICUs and hospitals.

This document is available to all distributors and also to nurses and NICUs.

Contact Us for More Information (800)676-0462

