Bebeonkers

Bebeonkers II (16-2) Our original patented design is an Oral/Nasal suction device that does not grab tissue and effectively removes the most tenacious secretions, quickly and without compromising the infant.

Now available: 16-0 Webeonkers for the ultra preemie infant, just like the 16-2 but smaller.

Packaged six per bag.





Flexible & Soft

Bebeonkers I (16-1) was designed for the Ultra Preemie infant to facilitate quick removal of secretions both nasally and orally but still retains a three-hole configuration for safe oral care without hurting baby.

Remember "Do No Harm"

Order #	Product Description	Bx Qty	
Bebeonkers Oral/Nasal Suction Device Only			
16-1	Bebeonkers I Ultra Preemie Bulbed Tip 6 pcs/zip bag,50 bgs/bx	300	
16-2	Bebeonkers II Standard Tip 6 pcs/zip bag, 50 bgs/bx	300	
16-0	Webeonkers Ultra Preemie STD tip 6 pcs/zip bag. 50 bags/bx	300	
16-1KK	Bebeonkers I with kover 3 Bebeonkers + 3 cover/bag, 50 bgs/bx	300	
16-2KK	Bebeonkers II with kover 3 Bebeonkers + 3 cover/bag, 50 bgs/bx	300	
16-0KK	Webeonkers with kover 3 Bebeonkers + 3 cover/bag, 50 bgs/bx	300	



Bebeonkers





Devices altered to facilitate the removal of oral secretions or designed ineffectively for tenacious secretions.

These devices were removed from hospitals around the country as Nurses make whatever they can to suction more effectively.

Major liability

Major liability

All above are single hole at distal end.

* remember a single hole device will grab tissue when suctioning and that's bad!!

Because there has never been an appropriate device for suctioning intubated or newborn infants, doctors, nurses and therapists have been using and creating inferior one-hole devices. They endanger the infants and make the hospital and medical professional liable if an infant is injured.


The patented
"BEBEONKERS" is that right device!!
It is the right tool for the right job.
Two or more holes and a tapered
inner wall does the job quickly and
without tearing tissue.



"Sometimes you just can't stand by when things are being done that can hurt a baby. We saw a problem that was international and made a product to solve it!" —Philip N. Rogone PA. RCP.

Bebe Keeper Kover


Can be used up to 12 hours on the same patient, or per hospital policy.





Bebe Keeper Kover is an effective cover to keep the Bebeonkers clean.

Can be used on Bebeonkers I and Bebeonkers II with suction on or off, protecting the device while not in use.



Order #	Product Description	Bx Qty		
Bebe Keeper Kover Only				
16-3	Bebe Keeper Kover for Bebeonkers 6 pcs/bag, 50 bgs/cs	300		

Can also be used with the Bebe-Ball.



Note new blue translucent color